

MesenCult™ Osteogenic Differentiation Kit (Human)



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #05465

1 Kit

For the *in vitro* differentiation of human MSCs into osteoblasts

Product Description

MesenCult™ Osteogenic Differentiation Kit (Human) is specifically formulated for the *in vitro* differentiation of human mesenchymal stromal cells (MSCs; also known as mesenchymal stem cells) into cells of the osteogenic lineage. This kit is suitable for the differentiation of human bone marrow (BM)- or adipose tissue (AT)-derived MSCs previously culture-expanded in serum-containing medium (e.g. MesenCult™ Proliferation Kit [Human; Catalog #05411] or MesenCult™-hPL Medium [Human; Catalog #05439]) or serum- and animal component-free MesenCult™-ACF Plus Medium (Catalog #05445). This kit is also suitable for the osteogenic differentiation of human embryonic stem (ES) and induced pluripotent stem (iPS) cell-derived mesenchymal progenitor cells (MPCs) generated with STEMdiff™ Mesenchymal Progenitor Kit (Catalog #05240).

NOTE: Complete MesenCult™ Osteogenic Differentiation Medium must be supplemented with L-Glutamine (Catalog #07100).

Product Information

The following components are sold as a complete kit (Catalog #05465) and are not available for individual sale.

COMPONENT NAME	COMPONENT #	SIZE	STORAGE	SHELF LIFE
MesenCult™ Osteogenic Differentiation Basal Medium (Human)	05466	200 mL	Store at 2 - 8°C.	Stable until expiry date (EXP) on label.
MesenCult™ Osteogenic Differentiation 5X Supplement (Human)	05467	50 mL	Store at -20°C.	Stable until expiry date (EXP) on label.

Preparation of Medium

Use sterile technique to prepare complete Osteogenic Differentiation Medium (Basal Medium + 5X Supplement + L-Glutamine). The following example is for preparing 50 mL of complete medium. If preparing other volumes, adjust accordingly.

NOTE: For aliquoting and storing either the supplement or the complete medium, polypropylene tubes are strongly recommended (e.g. Falcon® Conical Tubes, 15 mL [Catalog #38009]).

1. Thaw the 5X Supplement at room temperature (15 - 25°C). Alternatively, thaw supplement at 2 - 8°C overnight. Mix thoroughly.

NOTE: Once thawed, use immediately or aliquot and store at -20°C. Do not exceed the shelf life of the supplement. After thawing the aliquoted supplement, use immediately. Do not re-freeze.

2. Add 10 mL of 5X Supplement to 40 mL of Basal Medium. Mix thoroughly.
3. Add 0.5 mL of L-Glutamine (200 mM; Catalog #07100) to reach a final concentration of 2 mM. Mix thoroughly.

NOTE: If not used immediately, store complete Osteogenic Differentiation Medium at 2 - 8°C for up to 1 week. This medium does not contain antibiotics.

Directions for Use

Please read the entire protocol before proceeding.

For instructions on culturing human MSCs or ES/iPS-derived MPCs using the MesenCult™ media listed below, or instructions on derivation and expansion of cells with MPC-like properties from human ES or iPS cells using STEMdiff™ Mesenchymal Progenitor Kit (Catalog #05240), refer to the following Product Information Sheets available at www.stemcell.com.

- MesenCult™ Medium (Human; Catalog #05411)
- MesenCult™-hPL Medium (Human; Catalog #05439)
- MesenCult™-ACF Plus Medium (Catalog #05445)
- MesenCult™-ACF Plus Culture Kit (Catalog #05448)

For differentiating to the osteogenic lineage, use culture-expanded human MSCs between passages 1 - 4 and culture-expanded PSC-derived MPCs between passages 4 - 14.

The following protocol is for setting up differentiation assays using human BM- or AT-derived MSCs in a 6-well plate. If using other cultureware, adjust volumes accordingly.

NOTE: Only use tissue culture-treated cultureware.

1. Plate cells in 2 mL of growth medium per well. Recommended cell plating densities are as follows:

- MesenCult™ Medium (Catalog #05411): $8 - 16 \times 10^3$ cells/cm²
- MesenCult™-ACF Plus Medium or MesenCult™-hPL Medium (Human): $4 - 8 \times 10^3$ cells/cm²

NOTE: If using MesenCult™-ACF Plus Medium, ensure to coat cultureware as described in the Product Information Sheet (Document #10000003462).

2. Incubate at 37°C in 5% CO₂ until cells are approximately 80 - 90% confluent. This takes approximately 1 - 5 days depending on the expansion medium used.
3. Aspirate medium and replace with 2 mL of complete Osteogenic Differentiation Medium per well. Incubate at 37°C in 5% CO₂.
4. Change medium every 3 - 4 days until bone matrix formation occurs. The culture time for inducing differentiation is dependent on the cell source. See Table 1 for typical ranges.

Table 1. Culture Times for Inducing Osteogenic Differentiation for Various MSC Sources

MSC SOURCE	CULTURE TIMES FOR DIFFERENTIATION (DAYS)
Bone Marrow	10 - 15
Adipose Tissue	10 - 15
PSCs	14 - 21

5. Visualize osteogenic differentiation by staining with Alizarin Red S, alkaline phosphatase, or silver nitrate (von Kossa).

NOTE: The level of osteogenic differentiation for MSCs may vary depending on cell source, donor, and previous culture conditions.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2023 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, MesenCult, and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. Falcon is a registered trademark of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.