



EasyPlate™ EasySep™ Magnet

Magnet for Column-Free Immunomagnetic Separation

Processing up to 2×10^7 cells (mouse splenocytes) or
 1×10^7 cells (human PBMC) per well

Catalog # 18102



Scientists Helping Scientists™ | www.stemcell.com

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

The EasyPlate™ EasySep™ Magnet is an anodized aluminum body housing rare earth magnetic elements.

The EasyPlate™ EasySep™ Magnet is designed for cell separation procedures using selected EasySep™ reagents when processing multiple samples and a small number of cells (up to 2×10^7 cells per well). The EasyPlate™ Magnet generates a high-gradient magnetic field that is strong enough to separate cells labeled with magnetic EasySep™ particles without the use of columns. The EasyPlate™ EasySep™ Magnet is designed to hold a round-bottom non-tissue culture-treated 96-well plate.

Properties

Storage: Store at 15 - 25°C.

Shelf Life: Not applicable

STRONG MAGNETIC FIELD: Keep away from pacemakers, magnets, computer disks, watches, and other objects that respond to magnetic fields.

Handling / Directions For Use

The EasyPlate™ EasySep™ Magnet is only for use with certain EasySep™ kits. Follow the EasyPlate™ EasySep™ instructions provided with the EasySep™ cocktail of your choice. If your cell type or sample is not listed on www.stemcell.com, please contact techsupport@stemcell.com.

The EasyPlate™ EasySep™ Magnet is designed to hold a round-bottom non-tissue culture-treated 96-well plate (e.g. Costar Catalog #3788 or Corning Catalog #351177). If using a different type of non-tissue culture-treated 96-well plate, ensure that it fits properly on the EasyPlate™ EasySep™ magnet before use. Some 96-well plates may not sit flat on the magnet, which could affect the success of the separation.

CARE OF MAGNET: After completing a cell separation, wipe the magnet clean with distilled water and then with alcohol to dry the magnet. To disinfect the magnet, wipe with 10% diluted household bleach and then wipe off with a water-dampened cloth to avoid corrosion of metal parts. Do not submerge the magnet in liquid.



STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasyPlate, and EasySep are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.