

Carmine Alum



Dye for staining mouse mammary gland whole mounts

Catalog # 07070

100 mL

Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Carmine Alum is used and recommended for the whole mount staining of mouse mammary glands.

Properties

- Storage:** Store at 2 - 8°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Carmine and aluminum potassium sulfate in water
 - ProClin™ is added as a preservative

Handling / Directions For Use

1. Resect mammary glands from mice and spread onto a glass slide.
2. Fix in Carnoy's fixative for 2 - 4 hours at room temperature (15 - 25°C).
NOTE: Carnoy's fixative contains 60% ethanol, 30% chloroform, and 10% glacial acetic acid.
3. Wash in 70% ethanol for 15 minutes.
4. Change gradually to distilled water.
5. Wash once in distilled water for 5 minutes.
6. Change to Carmine Alum and stain overnight at room temperature (15 - 25°C).
7. Wash glands with sequential 15 minute washes in 70%, 95%, and 100% ethanol.
8. Clear with xylene and mount slides with Permount™ (Fisher Scientific Catalog #SP15-500).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. ProClin is a trademark of the Dow Chemical Company (Dow) or an affiliated company of Dow. Permount is a trademark of Fisher Scientific. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.