

Small Molecules

Sinefungin

Epigenetic modifier; Inhibits DNA methyltransferases (DNMTs)

Catalog # 73812

1 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713


INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Sinefungin is a DNA methyltransferase (DNMT) inhibitor ($IC_{50} = 0.1 - 20 \mu M$; Copeland et al.) and an analog of the methyl donor S-adenosyl-L-methionine (SAM; Horiuchi et al.; Nolan).

Molecular Name:	Sinefungin
Alternative Names:	A 9145; Antibiotic 32232RP; Antibiotic A 9145; RP 32232
CAS Number:	58944-73-3
Chemical Formula:	$C_{15}H_{23}N_7O_5$
Molecular Weight:	381.4 g/mol
Purity:	$\geq 95\%$
Chemical Name:	(2S,5S)-2,5-diamino-6-[(2R,3S,4R,5R)-5-(6-aminopurin-9-yl)-3,4-dihydroxyoxolan-2-yl]hexanoic acid
Structure:	



Properties

Physical Appearance:	A white solid
Storage:	Product stable at $-20^{\circ}C$ as supplied. Protect product from prolonged exposure to light. For long-term storage store with a desiccant. Stable as supplied for 12 months from date of receipt.
Solubility:	· Water ≤ 50 mM For example, to prepare a 10 mM stock solution in water, resuspend 1 mg in 262 μL of water.

Prepare stock solution fresh before use. Information regarding stability of small molecules in solution has rarely been reported, however, as a general guide we recommend storage in water at $-20^{\circ}C$. Aliquot into working volumes to avoid repeated freeze-thaw cycles. The effect of storage of stock solution on compound performance should be tested for each application.

For use as a cell culture supplement, stock solution should be diluted into culture medium immediately before use.

Published Applications

IMMUNOLOGY

- Antiparasitic activity against malarial, trypanosomal, and leishmanial species (Dube et al.; Nolan; Trager et al.).
- Inhibits Newcastle disease and vaccinia virus multiplication (Pugh et al.).
- Inhibits pneumococcal biofilm growth in vitro and colonization in an in vivo rat model (Yadav et al.).

References

- Copeland RA et al. (2009) Protein methyltransferases as a target class for drug discovery. *Nat Rev Drug Discov* 8(9): 724–32.
- Dube DK et al. (1983) Antitrypanosomal activity of sinefungin. *Am J Trop Med Hyg* 32(1): 31–3.
- Horiuchi KY et al. (2013) Assay development for histone methyltransferases. *Assay Drug Dev Technol* 11(4): 227–36.
- Nolan LL. (1987) Molecular target of the antileishmanial action of sinefungin. *Antimicrob Agents Chemother* 31(10): 1542–8.
- Pugh CS et al. (1978) Sinefungin, a potent inhibitor of virion mRNA(guanine-7-)-methyltransferase, mRNA(nucleoside-2'-)-methyltransferase, and viral multiplication. *J Biol Chem* 253(12): 4075–7.
- Trager W et al. (1980) *Plasmodium falciparum*: antimalarial activity in culture of sinefungin and other methylation inhibitors. *Exp Parasitol* 50(1): 83–9.
- Yadav MK et al. (2014) Sinefungin, a natural nucleoside analogue of S-adenosylmethionine, inhibits *Streptococcus pneumoniae* biofilm growth. *Biomed Res Int* 2014: 156987.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.