

ClonaCell™-HY Liquid HAT Selection Medium

Hybridoma selection medium, with HAT (serum-containing)

Catalog # 03831

500 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

ClonaCell™-HY Liquid HAT Selection Medium is a serum-containing medium with hypoxanthine, aminopterin, and thymidine (HAT) for the selection of hybridomas following the fusion of lymphocytes and myeloma cells. This medium has been verified for use in mouse and rat hybridoma development and reportedly is compatible for production and cloning of myelomas and/or hybridomas using lymphocytes from a variety of host animals including human, mouse, rat, and hamster.

- Simplifies hybridoma selection in liquid suspension cultures

Properties

- Storage:** Store at -20°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- DMEM
 - Serum
 - Hypoxanthine, aminopterin, and thymidine (HAT)
 - Gentamicin
 - 2-Mercaptoethanol
 - Phenol red
 - L-Glutamine and other supplements
 - Other ingredients

Handling / Directions For Use

1. Thaw ClonaCell™-HY Liquid HAT Selection Medium at room temperature (15 - 25°C) or overnight at 2-8°C. Mix well.
NOTE: Do not thaw ClonaCell™-HY Liquid HAT Selection Medium in a 37°C water bath.
2. If ClonaCell™-HY Liquid HAT Selection Medium is not used immediately, store at 2 - 8°C for up to 2 weeks. Alternatively, aliquot and store at -20°C until expiry date as indicated on the label.
3. Follow your own procedure for hybridoma generation, selection, and expansion or refer to the Technical Manual: ClonaCell™-HY: A Complete Workflow for Hybridoma Generation (available at www.stemcell.com or contact us to request a copy) for protocols using the complete range of ClonaCell™-HY products.

References

- Dorrell C et al. (2011) Isolation of mouse pancreatic alpha, beta, duct and acinar populations with cell surface markers. *Mol Cell Endocrinol* 339(1-2): 144-50.
- Fang L et al. (2008) Essential role of TNF receptor superfamily 25 (TNFRSF25) in the development of allergic lung inflammation. *J Exp Med* 205(5): 1037-48.
- García-Barreno B et al. (2014) Characterization of an enhanced antigenic change in the pandemic 2009 H1N1 influenza virus haemagglutinin. *J Gen Virol* 95(5): 1033-42.
- Jin C et al. (2005) Immunoglobulin G specifically binding plant N-glycans with high affinity could be generated in rabbits but not in mice. *Glycobiology* 16(4): 349-57.
- Weidanz JA et al. (2006) Levels of specific peptide-HLA class I complex predicts tumor cell susceptibility to CTL killing. *J Immunol* 177(8): 5088-97.
- Wittman VP et al. (2006) Antibody targeting to a class I MHC-peptide epitope promotes tumor cell death. *J Immunol* 177(6): 4187-95.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL QUALITY INFORMATION, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ClonaCell are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.