

EasySep™ Release Human PSC-Derived Neural Crest Cell Positive Selection Kit

For processing 1 x 10⁹ cells

Catalog #100-0047

Positive Selection

Document #10000007951 | Version 00



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Description

Isolate highly purified CD271+ cells from human pluripotent stem cell (PSC)-derived neural crest cell culture by immunomagnetic positive selection.

- Highly purified human PSC-derived neural crest cells isolated in less than 20 minutes
- No-wash removal of EasySep™ Releasable RapidSpheres™
- Compatible with downstream differentiation protocols

This kit targets human PSC-derived neural crest cells with antibody complexes recognizing the CD271 surface marker and EasySep™ Releasable RapidSpheres™. Desired cells are labeled with antibodies and magnetic particles, and separated without columns using an EasySep™ magnet. Unwanted cells are poured off, while desired cells remain in the tube. Then, bound magnetic particles are removed from the EasySep™-isolated human PSC-derived neural crest cells, which are immediately available for downstream applications such as flow cytometry or cell culture.

Following cell isolation with this EasySep™ Release kit, antibody complexes remain bound to the cell surface and may interact with Brilliant Violet™ antibody conjugates, polyethylene glycol-modified proteins, or other chemically related ligands.

Component Descriptions

COMPONENT NAME	COMPONENT #	QUANTITY	STORAGE	SHELF LIFE	FORMAT
EasySep™ Release Human PSC-Derived Neural Crest Cell Positive Selection Cocktail	300-0045	2 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A combination of monoclonal antibodies in PBS with 0.1% rHA.
EasySep™ Isolation Cocktail Enhancer	17900	1 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A solution that enhances the performance of the isolation cocktail.
EasySep™ Releasable RapidSpheres™ 50201	50201	2 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A suspension of magnetic particles in water.
EasySep™ Release Buffer (Concentrate)	20165	3 x 1 mL	Store at 2 - 8°C. Do not freeze.	Stable until expiry date (EXP) on label.	A buffer for release of Releasable RapidSpheres™ from cells following positive selection.

PBS - phosphate-buffered saline; rHA - recombinant human albumin

Components may be shipped at room temperature (15 - 25°C) but should be stored as indicated above.

Sample Preparation

Trypsin-EDTA (0.25%; Catalog #07901) is recommended for preparation of a single-cell suspension from human PSC-derived neural crest cells, as follows:

1. Discard culture medium. Rinse each well of a 6-well plate with 2 mL DMEM/F-12.
2. Add 1 mL Trypsin-EDTA (0.25%) to each well.
3. Incubate at 37°C for 5 - 7 minutes, until the cells are easily dislodged from the bottom of the plate with gentle tapping.
4. Add 1 mL Soybean Trypsin Inhibitor, ACF (Catalog #07457) to each well.
5. Pipette up and down 3 - 4 times with a 1 mL pipette tip to dissociate to single cells.
6. Combine the single-cell suspension from each well into a 15 mL conical tube (e.g. Catalog #38009).
7. Centrifuge at 300 x g for 10 minutes, with the brake on low.
8. Remove supernatant and resuspend cells at 2.5 x 10⁷ cells/mL in recommended medium.

NOTE: ACCUTASE™ (Catalog #07920) may be used in place of Trypsin. Refer to the Product Information Sheet for ACCUTASE™ for a cell detachment protocol.

Recommended Medium

RoboSep™ Buffer 2 (Catalog #20164) or PBS containing 0.5% bovine serum albumin (BSA) and 2 mM EDTA. Medium should be free of Ca⁺⁺ and Mg⁺⁺.

Directions for Use – Manual EasySep™ Protocol

See page 1 for Sample Preparation and Recommended Medium. Refer to Table 1 for detailed instructions regarding the EasySep™ procedure.

Table 1. EasySep™ Release Human PSC-Derived Neural Crest Cell Positive Selection Kit

STEP	INSTRUCTIONS	 EasySep™ (Catalog #18000)
1	Dilute Release Buffer (Concentrate) to prepare release buffer (1X).	Dilute 1 in 40 with recommended medium. NOTE: Release buffer (1X) must be prepared on the day of use. Refer to step 14 for required volume.
2	Prepare sample at the indicated cell concentration within the volume range.	2.5 x 10 ⁷ cells/mL 0.1 - 2 mL
	Add sample to required tube.	5 mL (12 x 75 mm) polystyrene round-bottom tube (e.g. Catalog #38007)
3	Add Selection Cocktail to sample.	50 µL/mL of sample
	Mix and incubate.	RT for 3 minutes
4	Vortex RapidSpheres™. NOTE: Particles should appear evenly dispersed.	30 seconds
5	Add Releasable RapidSpheres™ to sample.	50 µL/mL of sample
	Mix and incubate.	RT for 3 minutes
6	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL
	Place the tube (without lid) into the magnet and incubate.	RT for 1 minute
7	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring off the supernatant. Remove the tube from the magnet; this tube contains the isolated cells.	Discard supernatant
8	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL
9	Add Cocktail Enhancer to sample.	10 µL
	Mix and incubate.	RT for 3 minutes
10	Place the tube (without lid) into the magnet and incubate.	RT for 1 minute
11	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring off the supernatant. Remove the tube from the magnet; this tube contains the isolated cells.	Discard supernatant
12	Add recommended medium to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL
	Place the tube (without lid) into the magnet and incubate.	RT for 1 minute
13	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring off the supernatant. Remove the tube from the magnet; this tube contains the isolated cells.	Discard supernatant
14	Add release buffer (1X) to top up the sample to the indicated volume. Mix by gently pipetting up and down 2 - 3 times.	Top up to 2.5 mL
	Mix and incubate.	RT for 3 minutes
15	Place the tube (without lid) into the magnet and incubate.	RT for 1 minute
16	Pick up the magnet, and in one continuous motion invert the magnet and tube,* pouring the enriched cell suspension into a new tube.	Isolated cells (in the new tube) are ready for use

RT - room temperature (15 - 25°C)

* Leave the magnet and tube inverted for 2 - 3 seconds, then return upright. Do not shake or blot off any drops that may remain hanging from the mouth of the tube.

Notes and Tips

EASYSEPTM RELEASE BUFFER

EasySep™ Release Buffer (Concentrate) is supplied as a 40X concentrate; release buffer (1X) must be prepared on the day of use. To prepare release buffer (1X), dilute an appropriate volume 1 in 40 with recommended medium. Refer to step 14 of Table 1 for required volume.

ASSESSING PURITY

For purity assessment of human PSC-derived neural crest cells by flow cytometry, use the following fluorochrome-conjugated antibody clone:


- Anti-human CD271 antibody, clone ME20.4

NOTE: Viability is measured by exclusion of DAPI (Hydrochloride; Catalog #75004), Propidium Iodide (Catalog #75002), or 7-AAD (7-Aminoactinomycin D; Catalog #75001).


NOTE: Brilliant Violet™ antibody conjugates should be carefully titrated on EasySep™ Release-isolated cells prior to analysis by flow cytometry or fluorescence microscopy. For purity assessment with Brilliant Violet™ antibody conjugates, use of BD Horizon Brilliant™ Stain Buffer is recommended to reduce non-specific interactions. For more information, refer to the manufacturer's instructions or contact us at techsupport@stemcell.com.

Data


A H9 cells



B H7 cells



C 1C cells



Starting with a single-cell suspension of 1C, H7, or H9 cells differentiated using STEMdiff™ Neural Crest Differentiation Kit (Catalog #08610), the CD271+ cell content of the isolated fractions is typically $92.9 \pm 6.3\%$ (mean \pm SD using the purple EasySep™ Magnet).

(A) Starting with H9 cells, the purities of the start and final isolated fractions are 74.9% and 93.8%, respectively.

(B) Starting with H7 cells, the purities of the start and final isolated fractions are 47.3% and 90.8%, respectively.

(C) Starting with 1C cells, the purities of the start and final isolated fractions are 18.9% and 86.7%, respectively.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED. FOR ADDITIONAL INFORMATION ON QUALITY AT STEMCELL, REFER TO WWW.STEMCELL.COM/COMPLIANCE.

Copyright © 2020 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, RapidSpheres, and STEMdiff are trademarks of STEMCELL Technologies Canada Inc. ACCUTASE is a trademark of Innovative Cell Technologies. Brilliant Violet is a trademark of Sirigen Group Ltd. BD Horizon Brilliant is a trademark of Becton, Dickinson, and Company. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.