

SARS-CoV-2 (NS6) Peptide Pool

SARS-CoV-2 (NS6) peptide pool for immune cell activation

Catalog #100-0660

~25 µg/peptide



Scientists Helping Scientists™ | www.stemcell.com

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

The SARS-CoV-2 (NS6) Peptide Pool is a lyophilized mixture of 13 peptides from non-structural protein 6 (NS6) of SARS-CoV-2. SARS-CoV NS6 has been shown to localize in the endoplasmic reticulum/Golgi membrane of the host cell and disrupt nuclear import complex formation (Frieman et al.). This leads to the loss of STAT1 transport into the nucleus and inactivation of STAT1-activated genes that are required for antiviral effects (Bharath et al.; Frieman et al.). The sequence similarity of NS6 from the two coronaviruses suggests NS6 plays a similar role in SARS-CoV-2. The pool consists of 15-mer peptides with 11-amino-acid overlaps that cover amino acids 1 - 61 on NS6.

Product Information

Number of Peptides:	13
Source:	SARS-CoV-2 (severe acute respiratory syndrome coronavirus 2)
Protein ID:	P0DTC6 (Swiss-Prot)
Protein Name:	ORF6 Protein; Non-structural protein 6; NS6
Protein Sequence:	MFHLVDFQVTIAEILLIIMRTFKVSIWNLDYIINLIKNLSKSLTENKYSQLDEEQPMEID
Gene Name:	NS6
Purity:	Average 70%
Formulation:	Lyophilized as trifluoroacetate salts

Preparation and Storage

Storage:	Store at -20°C.
Stability:	Stable as supplied until expiry date (EXP) on label.
Preparation:	Warm to room temperature (15 - 25°C) before reconstitution. Add pure dimethyl sulfoxide (DMSO; ~40 µL) and dilute with water to the desired concentration. Final concentration of DMSO must be below 1% (v/v) to avoid toxicity in the biological system. If not used immediately, aliquot and store at -20°C. Protect from light. After thawing aliquots, do not re-freeze.

Related Products

For a complete list of peptide pools, as well as related products available from STEMCELL Technologies, visit www.stemcell.com or contact us at techsupport@stemcell.com.

References

Bharath BR et al. (2020) In silico screening of known small molecules to bind ACE2 specific RBD on spike glycoprotein of SARS-CoV-2 for repurposing against COVID-19. *F1000Res* 9: 663.

Frieman M et al. (2007) Severe acute respiratory syndrome coronavirus ORF6 antagonizes STAT1 function by sequestering nuclear import factors on the rough endoplasmic reticulum/Golgi membrane. *J Virol* 81(18): 9812–24.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.