

# Cytokines

## Mouse Recombinant RANKL, ACF

Receptor activator of nuclear factor kappa-B ligand, animal component-free



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog # 100-0943  
100-0944

100 µg  
1000 µg

## Product Description

Receptor activator of nuclear factor kappa-B ligand (RANKL) is a member of the tumor necrosis factor (TNF) superfamily (Anderson et al.). Cytokines in the TNF superfamily are involved in a variety of long-term cellular activities such as differentiation, proliferation, and cell death (MacEwan). RANKL is a type II homotrimeric transmembrane protein expressed in both a membrane-bound and secreted form (Ikeda et al.). RANKL binds to the receptor activator of nuclear factor kappa-B (RANK). Upon binding to its receptor, RANKL activates the AKT signaling pathway (Moon et al.). Osteoprotegerin (OPG) may also bind RANKL, and this binding competes with RANKL-RANK binding (Lacey et al.). RANKL is involved in osteoclastogenesis (Lacey et al.; Yasuda et al.) and T cell activation (Wong et al.). This product is animal component-free.

## Product Information

Alternative Names:	CD254, hRANKL2, ODF, OPGL, OPTB2, Osteoclast differentiation factor, soluble Receptor activator of NF-κB ligand, sOdf, TNF-related activation-induced cytokine, TNFSF11, TNF superfamily member 11, TNLG6B, Tumor necrosis factor superfamily member 11, TRAN
Accession Number:	O35235 (Pro143-Asp316) was expressed with an additional Met
Amino Acid Sequence:	MPAMMEGSWLD VAQRGKPEAQ PFAHLTINAA SIPSGSHKVT LSSWYHDRGW AKISNMTLSN GKLRVNQDGF YYLYANICFR HHETSGSVPT DYQLMVYV KTSIKIPSSH NLMKGGSTKN WSGNSEFHFY SINVGFFFKL RAGEEISIQV SNPSLLDPDQ DATYFGAFKV QDID
Predicted Molecular Mass:	19.5 kDa
Species:	Mouse
Formulation:	Lyophilized from a sterile-filtered aqueous solution containing 10 mM sodium phosphate and 50 mM sodium chloride, pH 7.5.
Source:	<i>E. coli</i>


## Specifications

Activity:	The EC50 of Mouse Recombinant RANKL, ACF is ≤ 50 ng/mL as determined by activity of RAW-Blue™ cells.
Purity:	≥ 95%
Endotoxin Level:	Measured by kinetic Limulus amoebocyte lysate (LAL) analysis and is ≤ 1.0 EU/µg protein.

## Preparation and Storage

Storage:	Store at -20°C to -80°C.
Stability:	Stable as supplied for 12 months from date of receipt.
Preparation:	Centrifuge vial before opening. Reconstitute the product in sterile water to at least 0.1 mg/mL by pipetting the solution down the sides of the vial. Do not vortex. The effect of storage of stock solution on product performance should be tested for each application. As a general guide, do not store at 2 - 8°C for more than 1 month or at -80°C for more than 3 months. Avoid repeated freeze-thaw cycles.

## Data



(A) The activity of Mouse Recombinant RANKL, ACF was tested by its ability to activate RAW-Blue™ cells. The EC<sub>50</sub> is defined as the effective concentration of the growth factor at which binding is at 50% of maximum. The EC<sub>50</sub> in the above example is 11.3 ng/mL.

(B) Mouse Recombinant RANKL, ACF was resolved with SDS-PAGE under reducing (+) and non-reducing (-) conditions and visualized by Coomassie Blue staining. Mouse Recombinant RANKL, ACF has a predicted molecular mass of 19.5 kDa (monomer).

## Related Products

For a complete list of cytokines, as well as related products available from STEMCELL Technologies, visit [www.stemcell.com/cytokines](http://www.stemcell.com/cytokines) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## References

- Anderson DM et al. (1997) A homologue of the TNF receptor and its ligand enhance T-cell growth and dendritic-cell function. *Nature* 390(6656): 175–9.
- Ikeda T et al. (2001) Determination of three isoforms of the receptor activator of nuclear factor-kappaB ligand and their differential expression in bone and thymus. *Endocrinology* 142(4): 1419–26.
- Lacey DL et al. (1998) Osteoprotegerin ligand is a cytokine that regulates osteoclast differentiation and activation. *Cell* 93(2): 165–76.
- MacEwan DJ. (2002) TNF ligands and receptors—a matter of life and death. *Br J Pharmacol* 135(4): 855–75.
- Moon JB et al. (2012) Akt induces osteoclast differentiation through regulating the GSK3β/NFATc1 signaling cascade. *J Immunol* 188(1): 163–9.
- Wong BR et al. (1997) TRANCE is a novel ligand of the tumor necrosis factor receptor family that activates c-Jun N-terminal kinase in T cells. *J Biol Chem* 272(40): 25190–4.
- Yasuda H et al. (1998) Osteoclast differentiation factor is a ligand for osteoprotegerin/osteoclastogenesis-inhibitory factor and is identical to TRANCE/RANKL. *Proc Natl Acad Sci USA* 95(7): 3597–602.

PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2021 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. RAW-Blue is a trademark of InvivoGen Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.