

Anti-Mouse Ly-6G Antibody, Clone 1A8

Antibodies

Rat monoclonal IgG2a antibody
against mouse Ly-6G, unconjugated



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Catalog #60031

500 µg 0.5 mg/mL

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

The 1A8 antibody reacts with Ly-6G, a 21 - 25 kDa GPI-anchored protein, which together with the structurally related Ly-6C protein comprises the granulocyte receptor-1 antigen (Gr-1). Gr-1 is expressed on monocytes, neutrophils and subsets of macrophages, plasmacytoid dendritic cells and T cells. Monocytes in the bone marrow transiently express Gr-1 during development and the expression level is strongly correlated with granulocyte differentiation and maturation. In the periphery, Gr-1 is found predominantly on neutrophils and is a useful marker for these cells. The 1A8 antibody binds specifically to Ly-6G, whereas another commonly used antibody, clone RB6-8C5, binds to both Ly-6G and Ly-6C. It has been reported that the 1A8 antibody detects Ly-6G-expressing granulocytes in peripheral blood, whereas the RB6-8C5 antibody also binds to Ly-6C-expressing lymphocytes, monocytes and dendritic cells.

Target Antigen Name:	Ly-6G
Alternative Names:	Gr-1, GR1, Ly6G
Gene ID:	546644
Species Reactivity:	Mouse
Host Species:	Rat (LEW)
Clonality:	Monoclonal
Clone:	1A8
Isotype:	IgG2a, kappa
Immunogen:	Mouse Ly-6G-transfected EL-4J cell line
Conjugate:	Unconjugated

Applications


Verified:	FC
Reported:	FA, FC, IP, IHC
Special Applications:	This antibody clone has been verified for purity assessments of cells isolated with EasySep™ kits, including EasySep™ Mouse Neutrophil Enrichment Kit (Catalog #19762).

Abbreviations: CellSep: Cell separation; ChIP: Chromatin immunoprecipitation; FA: Functional assay; FC: Flow cytometry; ICC: Immunocytochemistry; IF: Immunofluorescence microscopy; IHC: Immunohistochemistry; IP: Immunoprecipitation; WB: Western blotting

Properties

Formulation:	Phosphate-buffered solution, pH 7.2, containing 0.09% sodium azide
Purification:	The antibody was purified by affinity chromatography.
Stability and Storage:	Product stable at 2 - 8°C when stored undiluted. Do not freeze. For product expiry date, please request a lot-specific Certificate of Analysis from techsupport@stemcell.com .
Directions for Use:	For flow cytometry the suggested use of this antibody is ≤ 0.25 µg per 1 × 10 ⁶ cells in 100 µL volume. It is recommended that the antibody be titrated for optimal performance for each application.

Data



(A) Flow cytometry analysis of C57BL/6 mouse bone marrow cells labeled with Anti-Mouse Ly-6G Antibody, Clone 1A8, followed by anti-rat IgG, PE (black line histogram) or a rat IgG2a, kappa isotype control antibody followed by anti-rat IgG, PE (grey line histogram). Myeloid cells were gated for analysis.

(B) Immunohistochemical analysis of formalin-fixed, paraffin-embedded mouse uterine tissue labeled with Anti-Mouse Ly-6G Antibody, Clone 1A8, followed by anti-mouse, HRP. DAB substrate was used for visualization.

Related Products

PRODUCT NAME	CATALOG #	SIZE
Anti-Mouse Ly-6G Antibody, Clone 1A8	60031	500 μ g
Anti-Mouse Ly-6G Antibody, Clone 1A8, PE	60031PE	200 μ g
Anti-Mouse Ly-6G Antibody, Clone 1A8, PE	60031PE.1	50 μ g
Anti-Mouse Ly-6G Antibody, Clone 1A8, FITC	60031FI.1	50 μ g
Anti-Mouse Ly-6G Antibody, Clone 1A8, FITC	60031FI	500 μ g

References

1. Hestdal K, et al. Characterization and regulation of RB6-8C5 antigen expression on murine bone marrow cells. J Immunol 147(1): 22-28, 1991
2. Fleming TJ, et al. Selective expression of Ly-6G on myeloid lineage cells in mouse bone marrow. RB6-8C5 mAb to granulocyte-differentiation antigen (Gr-1) detects members of the Ly-6 family. J Immunol 151(5):2399-408, 1993 (IP)
3. Fleming TJ, Malek TR. Multiple glycosylphosphatidylinositol-anchored Ly-6 molecules and transmembrane Ly-6E mediate inhibition of IL-2 production. J Immunol 153(5): 1955-62, 1994 (FA)
4. Dietlin TA, et al. Mycobacteria-induced Gr-1+ subsets from distinct myeloid lineages have opposite effects on T cell expansion. J Leukocyte Biol 81(5): 1205-12, 2007 (FC)
5. Daley JM, et al. Use of Ly6G-specific monoclonal antibody to deplete neutrophils in mice. J. Leukocyte Biol 83(1): 64-70, 2008 (FC, Depletion)
6. Van Leeuwen, et al. Accumulation of myeloperoxidase-positive neutrophils in atherosclerotic lesions in LDLR^{-/-} mice. Arterioscler Thromb Vasc Biol 28(1): 84-89, 2008 (IHC)
7. Guiducci C, et al. Autoimmune skin inflammation is dependent on plasmacytoid dendritic cell activation by nucleic acids via TLR7 and TLR9. J Exp Med 207(13): 2931-42, 2010
8. Tadagavadi RK, Reeves WB. Endogenous IL-10 attenuates cisplatin nephrotoxicity: role of dendritic cells. J Immunol 185(8): 4904-11, 2010
9. Fujita M, et al. COX-2 blockade suppresses gliomagenesis by inhibiting myeloid-derived suppressor cells. Cancer Res 71(7): 2664-74, 2011
10. Rose S, et al. A novel Ly6C/Ly6G-based strategy to analyze the mouse splenic myeloid compartment. Cytometry A 81(4):343-50, 2012 (FC)

Copyright © 2013 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists and EasySep are trademarks of STEMCELL Technologies Inc. All other trademarks are the property of their respective holders.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.