

# CellAdhere™

## CellAdhere™ Dilution Buffer

Catalog #07183

100 mL



Scientists Helping Scientists™ | [www.stemcell.com](http://www.stemcell.com)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:info@stemcell.com) • [TECHSUPPORT@STEMCELL.COM](mailto:techsupport@stemcell.com)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

## Product Description

CellAdhere™ Dilution Buffer, when used in combination with Vitronectin XF™ (Catalog #07180), provides a defined surface for the culture of human embryonic stem (ES) cells and human induced pluripotent stem (iPS) cells.

## Properties

**Storage:** Store at 2 - 8°C.

**Shelf Life:** Stable until expiry date (EXP) on label.

**Contains:** CellAdhere™ Dilution Buffer is chemically defined.

## Directions for Use

Bring the CellAdhere™ Dilution Buffer to room temperature (15 - 25°C) prior to use.

For complete instructions on how to use CellAdhere™ Dilution Buffer with Vitronectin XF™, refer to the Vitronectin XF™ Product Information Sheet (Document #29285) available on our website at [www.stemcell.com](http://www.stemcell.com) or contact us to request a copy.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, CellAdhere, Vitronectin XF is developed and manufactured by Primorigen Biosciences Inc., and Vitronectin XF is a trademark of Primorigen Biosciences, Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.