

# MEM Non-Essential Amino Acid Solution (100X)

## MEM non-essential amino acid solution (100X)

Catalog # 07600

100 mL



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

MEM Non-Essential Amino Acid Solution (100X) has been pre-screened for use with other reagents of the ES-Cult™ product line for optimal performance for the in vitro hematopoietic differentiation of mouse embryonic stem (ES) cells and mouse induced pluripotent stem (iPS) cells.

## Properties

Storage:	Store at 2 - 8°C.
Shelf Life:	Stable until expiry date (EXP) on label.
Contains:	<ul style="list-style-type: none"><li>• 0.89 g/L L-alanine</li><li>• 1.5 g/L L-asparagine•H<sub>2</sub>O</li><li>• 1.33 g/L L-aspartic acid</li><li>• 1.47 g/L L-glutamic acid</li><li>• 0.75 g/L glycine</li><li>• 1.15 g/L L-proline</li><li>• 1.05 g/L L-serine</li></ul>

## Handling / Directions For Use

Use this product as directed in the protocol of choice.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2016 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, and ES-Cult are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.