

SpinSep™ Density Medium

Density Gradient Medium

Catalog # 17531

100 mL



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

FOR RESEARCH USE ONLY. NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES.

Product Description

SpinSep™ Density Medium is a density gradient medium. It is designed for the enrichment of specific human cell subsets.

Properties

- Storage:** Store at 15 - 25°C.
- Shelf Life:** Stable until expiry date (EXP) on label.
- Contains:**
- Iodixanol
 - Hetastarch
 - Sodium chloride
 - Sodium lactate
 - Dextrose
 - Calcium chloride
 - Potassium chloride
 - Magnesium chloride
 - Water

Handling / Directions For Use

Storage at 2 - 8°C is acceptable; ensure that the medium equilibrates to room temperature (15 - 25°C) and invert bottle to mix contents before use.

Use with the EasySep™ Human Myeloid Cell (CD33+ and CD66b+) Positive Selection Kit (Catalog #18653). Detailed information is included in the Product Information Sheet (Document #29009), available on our website at www.stemcell.com or contact us to request a copy.

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485 MEDICAL DEVICE STANDARDS.

Copyright © 2015 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, Scientists Helping Scientists, EasySep, and SpinSep are trademarks of STEMCELL Technologies Inc. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.