

Small Molecules

Rho Kinase Inhibitor IV

RHO/ROCK pathway inhibitor; Inhibits ROCK2

Catalog # 73802
73804

500 µg
1 mg



Scientists Helping Scientists™ | WWW.STEMCELL.COM

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

INFO@STEMCELL.COM • TECHSUPPORT@STEMCELL.COM


FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

Product Description

Rho Kinase Inhibitor IV is a selective and potent inhibitor of Rho-associated coiled-coil containing protein kinase 2 (ROCK2; $IC_{50} = 11.8$ nM; Tamura et al.). It is a glycyI analog of Fasudil (Catalog #73662) with increased specificity for ROCK2 (Tamura et al.). Rho Kinase Inhibitor IV is more potent than other ROCK inhibitors, including Y-27632 ($K_i = 220$ nM; Catalog #72302) and Fasudil ($IC_{50} = 158$ nM). It shows good specificity for ROCK2 compared to other kinases such as calcium/calmodulin-dependent kinase type II ($IC_{50} = 2.57$ µM), PKG ($IC_{50} = 2.35$ µM), Aurora A ($IC_{50} = 3.26$ µM), or PKA or PKC ($IC_{50} \geq 10$ µM each). ROCK1 and ROCK2 act downstream of the G protein Rho to regulate actin-myosin turnover and dynamics, and play an important role in stem cell renewal, smooth muscle contraction, cell adhesion, and proliferation (Narumiya et al.; Olson; Watanabe et al.). This product is supplied as the dihydrochloride salt of the molecule.

Molecular Name:	Rho Kinase Inhibitor IV (Dihydrochloride)
Alternative Names:	(S)-Glycyl-H-1152
CAS Number:	67-56-1
Chemical Formula:	$C_{18}H_{24}N_4O_3S \cdot 2HCl$
Molecular Weight:	449.4 g/mol
Purity:	$\geq 98\%$
Chemical Name:	2-amino-1-[(3S)-hexahydro-3-methyl-4-[(4-methyl-5-isoquinolinyl)sulfonyl]-1H-1,4-diazepin-1-yl]-ethanone, dihydrochloride

Structure:



Properties

Physical Appearance:	A solution in methanol
Storage:	Product stable at -20°C as supplied. Protect product from prolonged exposure to light. Stable as supplied for 12 months from date of receipt.
Solubility:	Not applicable.

Published Applications

DIFFERENTIATION

- Promotes neurite growth in primary rat neuronal cultures (Al-Ali et al.).
- Impairs primitive gut tube development including midgut elongation in *Xenopus* embryos (Reed et al.).

DISEASE MODELING

- Reduces intraocular pressure in a rabbit model of glaucoma (Tamura et al.).

References

- Al-Ali H et al. (2013) Chemical interrogation of the neuronal kinome using a primary cell-based screening assay. *ACS Chem Biol* 8(5): 1027–36.
- Narumiya S et al. (2009) Rho signaling, ROCK and mDia1, in transformation, metastasis and invasion. *Cancer Metastasis Rev* 28(1-2): 65–76.
- Olson MF. (2008) Applications for ROCK kinase inhibition. *Curr Opin Cell Biol* 20(2): 242–8.
- Reed RA et al. (2009) Morphogenesis of the primitive gut tube is generated by Rho/ROCK/myosin II-mediated endoderm rearrangements. *Dev Dyn* 238(12): 3111–25.
- Tamura M et al. (2005) Development of specific Rho-kinase inhibitors and their clinical application. *Biochim Biophys Acta* 1754(1-2): 245–52.
- Watanabe K et al. (2007) A ROCK inhibitor permits survival of dissociated human embryonic stem cells. *Nat Biotechnol* 25(6): 681–6.

Related Small Molecules

For a complete list of small molecules available from STEMCELL Technologies, visit www.stemcell.com/smallmolecules or contact us at techsupport@stemcell.com.

This product is hazardous. Please refer to the Safety Data Sheet (SDS).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.