

# Lambda® Plus Pipettor Starter Kit

Pipettor starter kit: 10 µL, 20 µL, 200 µL, 1000 µL

Catalog # 38063

1 Kit



Scientists Helping Scientists™ | [WWW.STEMCELL.COM](http://WWW.STEMCELL.COM)

TOLL FREE PHONE 1 800 667 0322 • PHONE +1 604 877 0713

[INFO@STEMCELL.COM](mailto:INFO@STEMCELL.COM) • [TECHSUPPORT@STEMCELL.COM](mailto:TECHSUPPORT@STEMCELL.COM)

FOR GLOBAL CONTACT DETAILS VISIT OUR WEBSITE

## Product Description

Corning® Lambda® Plus Pipettor is a single-channel, volumetric instrument designed to measure and transfer liquids precisely and safely. These fully autoclavable and UV-resistant pipettors have been designed for a universal fit with all common brands of pipette tips, including Corning® Filtered Pipette Tips (STEMCELL Catalog #38031). Each pipettor has smooth plunger movement, low pipetting force, and an adjustable tip ejector. Volume can be adjusted using the thumb wheel or plunger button, and a locking mechanism can be used to prevent accidental volume changes. The pipettor can accommodate volumes from 0.2 µL to 1000 µL, depending on the model. Pipettors are available individually and as a kit. The kit includes the 10 µL, 20 µL, 200 µL, and 1000 µL pipettors, as well as 10 µL, 200 µL, and 1000 µL pipette tips and a carousel stand.

2 µL Pipettor (Corning Catalog #4070) - Not included in Pipettor Starter Kit

10 µL Pipettor (Corning Catalog #4071)

20 µL Pipettor (Corning Catalog #4072)

200 µL Pipettor (Corning Catalog #4074)

1000 µL Pipettor (Corning Catalog #4075)

Pipettor Starter Kit (Corning Catalog #4069)

### APPLICATIONS

- Liquid handling
- RNA/DNA purification
- PCR
- Molecular biology

## Properties

**Storage:** Store at 15 - 25°C.

**Shelf Life:** Not applicable

## Features

### MATERIAL

- Lambda® Plus single-channel pipettors (one each of 10 µL, 20 µL, 200 µL, and 1000 µL formats)
- Universal carousel rack for six single-channel pipettors
- Pipette tips (one rack, 96 tips of each 10 µL, 200 µL, and 1000 µL formats)
- Calibration key
- Color-coded rings
- Tube of piston grease
- Identification stickers

## Handling / Directions For Use

For complete instructions, refer to the Technical Manual: Corning® Lambda® Plus Single-channel, 8-channel, and 12-channel Pipettors Instruction Manual (Document #DX21630), available at [www.stemcell.com](http://www.stemcell.com) or contact us to request a copy.

### WARRANTY

Three-year warranty is provided by Corning Life Sciences, valid from date of purchase. Original purchase receipt is required. For the complete warranty statement, refer to the Technical Manual or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Notes and Tips

For chemical resistance information, refer to the Technical Manual: Corning® Lambda® Plus Single-channel, 8-channel, and 12-channel Pipettors Instruction Manual or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

## Related Products

For a complete list of related products available from STEMCELL Technologies, visit [www.stemcell.com](http://www.stemcell.com) or contact us at [techsupport@stemcell.com](mailto:techsupport@stemcell.com).

STEMCELL TECHNOLOGIES INC.'S QUALITY MANAGEMENT SYSTEM IS CERTIFIED TO ISO 13485. PRODUCTS ARE FOR RESEARCH USE ONLY AND NOT INTENDED FOR HUMAN OR ANIMAL DIAGNOSTIC OR THERAPEUTIC USES UNLESS OTHERWISE STATED.

Copyright © 2018 by STEMCELL Technologies Inc. All rights reserved including graphics and images. STEMCELL Technologies & Design, STEMCELL Shield Design, and Scientists Helping Scientists are trademarks of STEMCELL Technologies Canada Inc. Corning and Lambda are registered trademarks of Corning Incorporated. All other trademarks are the property of their respective holders. While STEMCELL has made all reasonable efforts to ensure that the information provided by STEMCELL and its suppliers is correct, it makes no warranties or representations as to the accuracy or completeness of such information.